


## H<sub>2</sub>S Vial

### Always remember

1. To test the samples as soon as possible after collecting it from the source.
2. Mark the vial with sample ID and collection date with permanent marker before testing.

### Test Procedure:

1. Tare the wrap of the vial. Open the cap and fill the bottle upto the mark denoted on label.
2. Close the vial and invert a couple of time to mix the contents.
3. Incubate the vial at room temperature for 24-30 hours.
4. After incubation look for the colour change inside the vial from yellow to black.
5. Blackening of the sample or the paper inside denotes presence of H<sub>2</sub>S producing bacteria in the sample. No change in colour indicates absence of H<sub>2</sub>S producing bacteria.



**Store at 15-20°C in dry and dark place**